

Remarkable separations guaranteed time after time







Syncronis HPLC and UHPLC Columns

Remarkable separations guaranteed time after time

When developing a new method, one of the most important goals for the chromatographer is to achieve a consistent, reproducible separation. The selection of a highly reproducible HPLC column is essential if this goal is to be attained.


Our Thermo Scientific[™] Syncronis[™] HPLC columns are manufactured, packed and tested in ISO9000:2008 accredited facilities. Each lot of silica is tested for the physical properties of the silica support and only released for production if it meets the stringent test specifications.

Each bonded lot of chromatographic packing material is rigorously tested for primary and secondary interactions with the bonded phase.

New, enhanced, automated packing methods drive consistency even further and every column is individually tested to ensure that it meets the required quality.

These extensive testing and quality control procedures ensure the delivery of a consistent product, column after column.





Syncronis HPLC columns are available in a range of chemistries to give reproducible separations in reversed phase, HILIC and normal phase chromatography.

Testing and Quality Control

Silica Characterization	4
Bonded Phase Characterization	5
Column Packing	6
Summary of Tests Performed	6

Syncronis Stationary Phases

1.7 µm Particles for UHPLC Applications	7
C18	8–9
C8	10–11
aQ	12–13
Phenyl	14
Amino	15
Silica	16
HILIC	17
Ordering Information	18

Testing and Quality Control

Silica Characterization

Consistent separations require a rugged, reproducible silica backbone and to achieve this consistency, strict control of the physical properties of the silica particles is essential. Each manufactured lot of Syncronis silica is tightly controlled and extensively tested for its particle size and distribution, pore size and surface area and elemental purity. Only those batches which meet the rigorous quality control specifications are used.

Particle size and distribution





Tight control of the particle classification process ensures that a narrow particle size distribution is achieved around the target particle size, an important consideration for consistent chromatographic efficiency. Each manufactured lot of Syncronis silica is tested for its particle size and distribution using a laser particle size analyser. Three particle sizes are available: 1.7 μm for rapid UHPLC separations plus 3 μm and 5 μm for the more traditional HPLC analysis.

Pore size and surface area





Well controlled pore size and surface area are key to ensuring consistent carbon load and retentive properties of the chromatographic media. Each batch of Syncronis silica is tested for its pore size and surface area using liquid nitrogen adsorption.

Syncronis columns are based on highly pure 100 Å silica, with a surface area of 320 m^2/g , compared to 200 m^2/g for typical silica based material. This greater surface area ensures good retention of analytes having a range of hydrophobicity, away from the solvent front. The high surface area also allows for higher sample loading.

Silica purity (metals content)

The purity of the silica support is of particular importance when considering the separation of polar and basic compounds. Older, less pure silica supports contain a greater number of metallic impurities. The presence of certain metallic impurities with electron withdrawing properties (particularly aluminium) in silica can activate the silanols so that they become highly acidic, which can lead to peak tailing for basic solutes. Metallic impurities can also complex with chelating solutes, resulting in asymmetrical or tailing peaks. In extreme cases, these interactions may be strong enough to result in complete retention of the solute.

Each batch of Syncronis silica is tested for metals content using atomic emission spectroscopy.

Bonded Phase Characterization

Syncronis HPLC columns are bonded and endcapped with a range of stationary phases to effect different selectivity in separation. Whatever the bonded phase, rigorous testing and precise control of the bonding process are essential to achieve consistent chromatography. Each batch of Syncronis chromatographic media is tested for carbon load and characterized by stringent chromatographic testing before it is used to pack columns.

Carbon loading and surface coverage





The hydrophobic retention of a stationary phase is directly dependent on the carbon loading on the silica. Precise control of the batch to batch carbon loading is therefore a critical factor in ensuring consistent retention times. Syncronis reversed phase columns are densely bonded and double endcapped to minimise the number of residual silanols available to interact with basic analytes. Each batch of bonded phase is tested for carbon load using a total carbon analyzer.

Chromatographic tests

The retention properties of a reversed-phase packing material can be categorized into hydrophobic interactions, which include the measure of the hydrophobicity of the ligand and its density, steric or shape selectivity and secondary interactions such as silanol and surface metal activity. The impact that interactions between analytes and silanols have on the chromatographic performance depends on the pH of the mobile phase. Silanols on the silica surface can hydrogen bond (both as a donor and acceptor) and dissociated silanols can ion exchange with protonated bases.

To ensure consistent, predictable separations, the chromatographic media packed into Syncronis HPLC columns is extensively characterized using a series of diagnostic chromatographic tests, based on those developed by Tanaka¹. These tests rigorously probe interactions between analytes and stationary phase, measuring hydrophobicity, shape selectivity and secondary interactions with bases, acids and chelators.

1. K. Kimata, K. Iwaguchi, S. Onishi, K. Jinno, R. Eksteen, K. Hosoya, M. Araki and N. Tanaka, J.Chromatographic Science, v. 27 p. 721-728 (1989)

Column Packing

Every Syncronis HPLC column is individually tested and will not be released unless it meets the required retention, efficiency and peak symmetry. This testing is used to confirm the quality of the column packing process and the stability of the packed bed inside the column.

To ensure the most consistent results column after column, all Syncronis HPLC columns are packed using automated workstations.

To ensure that the measurements are a true indication of the quality of the packing, the testing of individual columns is performed on a highly optimized system.

Summary of Tests Performed on Syncronis HPLC Columns

	Test	C18	83	aQ	Phenyl	Amino	Silica	HILIC
Silica	Particle size and distribution	Y	Y	Y	Υ	Υ	Y	Υ
	Pore size and surface area	Y	Υ	Υ	Υ	Υ	Υ	Υ
Bonded phase	Carbon Load	Y	Υ	Υ	Υ	Υ		Υ
Chromatographic	Hydrophobic retention	Y	Υ	Υ	Υ			
	Hydrophobic selectivity	Υ	Υ	Υ	Υ			
	Steric selectivity	Υ	Υ	Υ	Υ			
	Hydrogen bonding capacity	Υ	Υ	Υ	Υ			
	Activity towards basic compounds	Υ	Υ	Υ	Υ			
	Ion exchange capacity (pH 7.6)	Y	Υ	Υ	Υ			
	Activity towards acidic compounds	Υ	Υ	Υ	Υ			
	Ion exchange capacity (pH 2.7)	Υ	Υ	Υ	Υ			
	Anion exchange test					Υ		
	Normal phase test						Υ	
	HILIC retention and selectivity							Υ
Column packing	Reversed phase packing test	Y	Υ	Υ	Υ			
	Normal phase packing test					Υ	Υ	
	HILIC packing test							Υ

1.7 µm Particles for UHPLC Applications

 $1.7 \mu m$ particles give higher efficiency than 3 μm or 5 μm particles and this efficiency is delivered over a greater range of optimum linear velocity. This makes it possible to operate at higher flow rates without losing performance. Because shorter columns packed with $1.7 \mu m$ particles give equivalent efficiency to longer columns packed with $5 \mu m$ particles faster analysis and solvent savings for the chromatographer become a reality.

Three Tips for Method Transfer

- To maintain an equivalent separation when transferring a method it is important to keep the linear velocity constant between the original and new method.
- 2. Sub-2 μm-based methods are most often transferred to smaller volume columns, so the same injection volume will take up a larger proportion of the new column, possibly leading to band broadening. It is therefore important to scale down the injection volume to match the change in column volume.
- 3. Geometrical transfer of the gradient requires calculation of the number of column volumes of mobile phase in each segment (time interval) of the gradient in the original method to ensure that the new calculated gradient takes place over the same number of column volumes, for the new column.

We also offer a convenient HPLC method transfer calculator at the Chromatography Resource Center www.thermofisher.com/SBE

System Considerations

With 1.7 µm particles, analyses can be performed with a high linear velocity through the column without loss in performance, provided the LC system is optimized to operate under these conditions. In order to produce fast, efficient chromatography, all system components for the assay should also be considered. Modern ultra high pressure liquid chromatography (UHPLC) instruments, including the Thermo Scientific Vanquish UHPLC System will take account of these factors.

There are three major system considerations to remember when using short columns packed with 1.7 µm particles.

- 1. The system volume (connecting tubing ID and length, injection volume, UV detector flow cell volume) must be minimized.
- The detector time constant and sampling rate need to be carefully selected.
- When running fast gradients pump delay volume needs to be minimal.




Syncronis C18


Outstanding column to column reproducibility

When developing a new method, one of the most important goals for the chromatographer is to achieve a consistent, reproducible separation. The selection of a highly reproducible HPLC column is essential if this goal is to be attained.

- · Outstanding reproducibility
- Highly pure, high surface area silica (320 m²/g)
- · High carbon load for increased retention
- · Double endcapped for extra surface coverage
- · Highly inert towards basic compounds
- · Rigorously tested to ensure quality



Syncronis C18 columns show excellent column to column reproducibility



Column: Syncronis C18 5 μ m, 250 mm \times 4.0 mm Mobile Phase: water:methanol (4:1)

Flow Rate: 1.0 mL/min Inj. Volume: 10 µL
Temp.: 25 °C
Detection: 265 nm


Sample: Zidovudine

Column	Retention Time (min)	Efficiency	Peak Area
1	25.82	62069	11532105
2	26.03	61688	11543904
3	25.90	62657	11527718
4	25.66	61317	11463444
5	25.92	63142	11520618
Inter-column precision (% RSD)	0.52	1.18	0.27

Syncronis C18 HPLC columns show excellent column to column reproducibility, as illustrated here by the analysis of zidovudine using five separate columns. The reproducibility in terms of retention time and peak area is less than or equal to 0.5%, column to column. This indicates that the columns are well packed.

The variation in peak area is 0.27%, which is important for quantitation of analytes.

Application: Chloramphenicol (USP)



Column:Mobile Phase: Flow Rate: Inj. Volume:

Inj. Volume: Temp.: Detection:

Sample:

Syncronis C18 5 μ m, 100 mm \times 4.6 mm


water:methanol:glacial acetic acid (54.9:45:0.1) 1.0 mL/min 10 μL 25 °C

Chloramphenicol

280 nm

Parameter	USP specification	Measured (6 replicate injections)		
Efficiency (N)	> 1800	6164		
Tailing factor	< 2	1.06		
%RSD retention time	< 1%	0.03%		
%RSD peak area	< 1%	0.32%		

Application: Ibuprofen and Valerophenone (USP)



Column:Mobile Phase: Flow Rate:

Inj. Volume: Temp.: Detection:

Sample:

Syncronis C18 5 μ m, 150 mm \times 4.0 mm

water/phosphoric acid (pH 2.5):acetonitrile (66.3:33.7) 2.0 mL/min $\,$

5 μL 30 °C 214 nm

1. Valerophenone

2. Ibuprofen


Parameter	USP Specification	Measured Valerophenone (5 replicate injections)	Measured Ibuprofen (5 replicate injections)
Resolution	> 2.0	_	3.56
Relative retention time	~ 0.8	0.84	_
Efficiency (N)	_	8317	5872
Tailing factor	_	1.11	1.38
%RSD retention time	_	0.50%	0.77%
%RSD peak area	_	0.50%	0.27%

Syncronis C8


Low hydrophobicity columns for less retention than Syncronis C18

Syncronis C8 HPLC columns are less hydrophobic than the C18 and are therefore particularly useful where the lesser degree of hydrophobicity is needed in order to successfully retain compounds of interest. Syncronis C8 HPLC columns can also be used where it is desirable to elute compounds more quickly.

- Outstanding reproducibility
- Highly pure, high surface area silica (320 m²/g)
- Less hydrophobic than C18
- · Double endcapped for extra surface coverage
- · Highly inert towards basic compounds
- · Rigorously tested to ensure quality



Application: Fenoprofen (USP)




Column: Syncronis C8 5 μm, 150 mm × 4.6 mm Mobile Phase: acetonitrile:water:phosphoric acid (50:49.6:0.4) Flow Rate: 2.0 mL/min Inj. Volume: 20 μL Temp.: 30 °C Detection: 272 nm

2. Gemfibrozil

Sample: 1. Fenoprofen

Parameter	USP Specification	Measured Fenoprofen (5 replicate injections)	Measured Gemfibrozil (5 replicate injections)
Resolution	> 8	-	17.6
Relative retention time	~ 0.5	0.48	_
Efficiency (N)	> 3000	9812	10254
Tailing factor	< 2	1.21	1.22
%RSD retention time	< 2%	0.13%	0.14%
%RSD peak area	< 2%	1.6%	1.8%

Application: Uron Herbicides



Column: Syncronis C8 5 μ m, 150 mm \times 4.6 mm

Mobile Phase:

A: water
B: acetonitrile
35 to 60 % B in 10 minutes
1.0 mL/min Gradient:

Flow Rate: 20 μL 30 °C 240 nm Inj. Volume: Temp.: Detection:

Sample: 1. Tebuthiuron

2. Metoxuron

3. Monuron

4. Chlorotoluron
5. Diuron

6. Linuron


Herbicide	RT (%RSD) (6 replicate injections)	Peak Area (%RSD) (6 replicate injections)	Peak Asymmetry
1 - Tebuthiuron	0.31	0.95	1.17
2 - Metoxuron	0.25	0.64	1.18
3 - Monuron	0.18	0.20	1.16
4 - Chlorotoluron	0.12	0.55	1.15
5 - Diuron	0.10	0.37	1.19
6 - Linuron	0.05	0.65	1.13

Syncronis aQ


Controlled interaction to retain and resolve polar analytes

In comparison to a conventionally endcapped C18, the Syncronis aQ polar end-capped C18 stationary phase exhibits superior stability towards aqueous mobile phase. Syncronis aQ shows no degradation in performance after 100 injections in a buffered 100 % aqueous eluent.

- Stability in 100 % aqueous mobile phase
- · Outstanding reproducibility
- Highly pure, high surface area silica (320 m²/g)
- · Highly inert towards basic compounds
- · Rigorously tested to ensure quality



Stability in 100% aqueous mobile phase



Column: Mobile Phase: Flow Rate: Inj. Volume: Temp.:

Detection:

Sample:

Syncronis aQ 5 µm, 100 mm × 4.6 mm

50 mM aqueous potassium hydrogen phosphate (pH 6) 0.7 ml /min


2 μL 30 °C 260 nm

1. Cytidine-5'-diphosphate


Adenosine-5'-triphosphate

3. Adenosine-5'-diphosphate

4. Adenosine-5'-monophosphate




In contrast, the performance of the C18 packing begins to deteriorate appreciably after roughly 50-60 replicate injections of the mixture of analytes. The decline in chromatographic performance is more pronounced for the later-eluting compounds. As shown on the right, there is a 20 % decrease in retention time for adenosine monophosphate on the C18 column.



Comparison of relative retention time for 5-AMP on Syncronis aQ and C18 over 100 injections


Application: Lamivudine (USP)



Column: Syncronis aQ 5 μ m, 250 mm × 4.6 mm

Sample: Lamivudine

Application: Amoxicillin and Potassium Clavulanate (USP)



Column:

Mobile Phase: Flow Rate: Inj. Volume: Temp.: Detection:

Sample:

210 nm

2.0 mL/min 20 µL 25 °C

Amoxicillin
 Potassium Clavulanate

Syncronis aQ 5 μ m, 300 mm imes 4.0 mm

phosphate buffer (pH 4.4):methanol (95:5)

Parameter	USP Specification	Measured Amoxicillin (6 replicate injections)	Measured K Clavulanate (6 replicate injections)
Resolution	> 3.5	_	12.8
Efficiency (N)	> 550	7598	6475
Tailing factor	< 1.5	1.15	0.92
%RSD retention time	< 2%	0.29%	0.36%
%RSD peak area	< 2%	0.30%	0.29%

Syncronis Phenyl


Enhanced retention of aromatic compounds

Syncronis Phenyl HPLC columns provide alternative selectivity to C18 and are particularly useful for the retention and separation of aromatic compounds.

- · Outstanding reproducibility
- Highly pure, high surface area silica (320 m²/g)
- Alternative selectivity to C18
- Double endcapped for extra surface coverage
- · Highly inert towards basic compounds
- · Rigorously tested to ensure quality

PH Range 0 2 to 8 14 Hydrophobicity Low High Pore Size 0 100Å 300 Carbon Load (%) 0 11 25 Particle Size 1.7μm, 3μm, 5μm USP L11

Application: Oxacillin Sodium (USP)



Parameter	USP Specification	Measured (6 replicate injections)
Efficiency (N)	> 2000	7904
Tailing factor	< 1.6	1.42
%RSD retention time	< 2%	0.03%
%RSD peak area	< 2%	0.29%

 $\begin{array}{ll} \text{Inj. Volume:} & 10 \ \mu\text{L} \\ \text{Temp.:} & 25 \ ^{\circ}\text{C} \\ \text{Detection:} & 225 \ \text{nm} \end{array}$

Sample: Oxacillin Sodium (0.11mg/mL)

Syncronis Amino


Versatile aminopropyl phase

Syncronis Amino HPLC columns exhibit excellent chromatographic properties in weak anion exchange, reversed phase, normal phase and HILIC.

- Outstanding reproducibility for reversed phase, normal phase, ion exhange and HILIC
- Highly pure, high surface area silica (320 m²/g)
- Alternative selectivity to C18
- Double endcapped for extra surface coverage
- · Rigorously tested to ensure quality

Per Range 0 2 to 8 14 Hydrophobicity Low High Pore Size 0 100Å 300 Carbon Load (%) 0 4 25 Particle Size 1.7μm, 3μm, 5μm USP L8

Application: Lactulose



Column: Syncronis Amino 5 μ m, 150 mm \times 4.6 mm

Mobile Phase: water:acetonitrile (30:70)

Flow Rate: 1.0 mL/min Inj. Volume: 5 μ L Temp.: 35 °C Detection: RI

Sample: Lactulose




Syncronis Silica


For highly efficient, normal phase chromatography

Syncronis Silica HPLC columns serve as a powerful and efficient tool for the chromatography of moderately polar organic compounds by normal phase chromatography.

- Highly pure, high surface area silica (320 m²/g)
- · Excellent reproducibility for normal phase chromatography
- · Rigorously tested to ensure quality



Application: Cetirizine (USP)



Parameter	USP Specification	Measured (6 replicate injections)	
Tailing factor	< 2.0	1.05	
%RSD peak area	< 2%	0.17%	

Column: Syncronis Silica 5 µm, 250 mm x 4.6 mm Mobile Phase: acetonitrile:water:sulfuric acid (93:6.6:0.4) Flow Rate: 1.0 mL/min

| Inj. Volume: 10 μL | Temp.: 30 °C | Detection: 230 nm


Sample: Cetirizine

Syncronis HILIC

Enhanced retention of polar and hydrophilic analytes


Syncronis HILIC is based on highly pure, high surface area silica particles. The zwitterionic modified stationary phase results in total charge equalisation and therefore a neutral (uncharged) but highly polar surface. Syncronis HILIC HPLC columns offer enhanced retention of polar and hydrophilic analytes.

- Highly pure, high surface area silica (320 m²/g)
- · Zwitterionic bonded phase
- · Enhanced retention of polar and hydrophilic analytes
- · Excellent reproducibility
- · Rapid equilibration
- · Rigorously tested to ensure quality



Particle Size 1.7µm, 3µm, 5µm

Application: Allantoin



Column: Syncronis HILIC 5 µm, 100 mm × 4.6 mm

Mobile Phase: ammonium formate buffer (pH 3):acetonitrile (10:90) Flow Rate: 1.0 mL/min Inj. Volume: 10 μ L Temp.: 30 °C

Detection: 210 nm

Sample: Allantoin

Ordering Information

Syncronis HPLC Columns

Particle Size (μm)	Length (mm)	ID (mm)	C18	C8	aQ	Phenyl	Amino	Silica	HILIC
1.7	30	2.1	97102-032130	-	_	-	97702-032130	_	-
	50	2.1	97102-052130	97202-052130	97302-052130	97902-052130	97702-052130	97002-052130	97502-052130
	50	3.0	97102-053030	_	_	_	_	_	_
	50	4.6	97102-054630	97202-054630	97302-054630	97902-054630	97702-054630	97002-054630	97502-054630
	100	2.1	97102-102130	97202-102130	97302-102130	97902-102130	97702-102130	97002-102130	97502-102130
	100	3.0	97102-103030	97202-103030	97302-103030	97902-103030	97702-103030	97002-103030	97502-103030
3	30	2.1	97103-032130	-	_	_	_	_	-
	50	2.1	97103-052130	97203-052130	97303-052130	97903-052130	97703-052130	97003-052130	97503-052130
	50	3.0	97103-053030	_	_	_	_	_	_
	50	4.6	97103-054630	-	_	-	-	_	-
	100	2.1	97103-102130	-	_	_	_	_	-
	100	3.0	97103-103030	97203-103030	97303-103030	_	_	_	_
	100	4.6	97103-104630	97203-104630	97303-104630	97903-104630	97703-104630	97003-104630	97503-104630
	150	2.1	97103-152130	-	_	-	_	_	-
	150	4.6	97103-154630	97203-154630	97303-154630	97903-154630	97703-154630	97003-154630	97503-154630
5	30	2.1	97105-032130	97205-032130	97305-032130	97905-032130	97705-032130	97005-032130	97505-032130
	50	2.1	97105-052130	97205-052130	97305-052130	97905-052130	97705-052130	97005-052130	97505-052130
	50	3.0	97105-053030	97205-053030	97305-053030	97905-053030	97705-053030	97005-053030	97505-053030
	50	4.6	97105-054630	97205-054630	97305-054630	97905-054630	97705-054630	97005-054630	97505-054630
	100	2.1	97105-102130	97205-102130	97305-102130	97905-102130	97705-102130	97005-102130	97505-102130
	100	3.0	97105-103030	97205-103030	97305-103030	97905-103030	97705-103030	97005-103030	97505-103030
	100	4.0	_	97205-104030	_	_	_	_	_
	100	4.6	97105-104630	97205-104630	97305-104630	97905-104630	97705-104630	97005-104630	97505-104630
	150	2.1	97105-152130	97205-152130	97305-152130	97905-152130	97705-152130	97005-152130	97505-152130
	150	3.0	97105-153030	97205-153030	97305-153030	97905-153030	97705-153030	97005-153030	97505-153030
	150	4.0	97105-154030	97205-154030	97305-154030	97905-154030	97705-154030	97005-154030	97505-154030
	150	4.6	97105-154630	97205-154630	97305-154630	97905-154630	97705-154630	97005-154630	97505-154630
	250	2.1	97105-252130	-	_	-	_	-	-
	250	3.0	97105-253030	97205-253030	97305-253030	97905-253030	97705-253030	97005-253030	97505-253030
	250	4.0	97105-254030	97205-254030	97305-254030	97905-254030	97705-254030	97005-254030	97505-254030
	250	4.6	97105-254630	97205-254630	97305-254630	97905-254630	97705-254630	97005-254630	97505-254630

Other column dimensions are available. Please contact Customer Services for more details

Syncronis HPLC Guard Columns (5 µm Particle Size)

Length (mm)	ID (mm)	Guard Holder	C18	C8	aQ	Phenyl	Amino	Silica	HILIC
10	2.1	852-00	97105-012101	97205-012101	97305-012101	97905-012101	97705-012101	97005-012101	97505-012101
10	3.0	852-00	97105-013001	97205-013001	97305-013001	97905-013001	97705-013001	97005-013001	97505-013001
10	4.0/4.6	850-00	97105-014001	97205-014001	97305-014001	97905-014001	97705-014001	97005-014001	97505-014001

Resources

for Chromatographers

Thermo Scientific Chromatography Columns and Consumables Catalog

This extensive catalog offers 450 pages of proven chromatography tools and product selection guides. Available online, with a robust search tool and optimized for your iPad®.

Visit www.thermofisher.com/catalog



Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit www.thermofisher.com/SBE

ChromeExpert

A dedicated team of technical support specialists and application chemists available to help you with product selection and assistance when using your chromatography consumables. Our experts have access to the latest chromatography technology and will act as your trusted advisors. For application and technical support.

Visit www.thermofisher.com/chromexpert

AppsLab Library

Thermo Scientific™ AppsLab Library of Analytical Applications provides more than 1300 detailed application examples for the columns listed in the 2016–2017 Chromatography Columns and Consumables Catalog. Search, filter and download complete methods to optimize your separation or implement validated methods using Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System (CDS) software. AppsLab Library makes our global application expertise accessible to you—online and downloadable.

Visit www.thermofisher.com/AppsLab









